

PDM module User Manual

1.0 Powering up the module

In order to switch on the module, simply plug in the AC adapter in the power supply connector. See *Table 1* for input DC Electrical requirements. After having applied the power, the module will switch on and, while all the internal diagnostic tests are performed, it will turn on the diagnostic output LEDs with green colour. After 20-25 seconds, if no error will occur, the LEDs will be switched off and the module will be ready for photon counting. If an error is encountered, all the internal voltages will be switched off and the diagnostic LEDs will display the problem using the colour code explained in *Table 2*. The module will remain in this situation until the power is removed.

Table 1. DC Input electrical characteristics

	<i>Min</i>	<i>Typ</i>	<i>Max</i>	<i>Unit</i>
<i>Input voltage</i>	5		12	Volt
<i>Input Power</i>			7	Watt

Note: input absolute maximum rating is 16V

Table 2. Diagnostic LED diagram

	TOP LED	BOTTOM LED
Vin higher than 13V	Red	Off
Vin lower than 4.5V	Orange	Off
Wrong SPAD Bias Voltage	Orange	Orange
SPAD parameters memory fault	Blinking Orange	Blinking Orange
SPAD working temperature not reached or not stable	Off	Orange
Peltier current saturated	Off	Red
Normal Condition	Green for 20-25 seconds after power-up, then Off	Green for 20-25 seconds after power-up, then Off

2.0 Inputs – Outputs

The module has one standard TTL output and one standard TTL gate input. Modules with the additional timing module installed have a second timing output.

TTL OUT

output: TTL OUT is designed to drive 50 Ω coaxial cables and is internally series terminated as conceptually shown in *Figure 1*. This output is short circuit protected. The TTL output connector is a standard BNC. TTL OUT low level is 0V and TTL OUT high level when driving a 50 Ω terminated input is about 3.5V.

Timing

output: NIM-standard fast-negative output. Typical NIM pulse amplitude is -800mV on a 50 Ω terminated transmission line. In order to obtain the best performances, it is recommended to use the NIM timing output only with terminated 50 Ω coax cables. NIM output connector is a LEMO 50 Ω NIM-CAMAC series (ERN.00.250.CTL).

GATE IN

Input: Gate-in input is designed to accept a standard TTL signal, AC terminated. The impedance is 10k Ω in DC and 50 Ω on the rising and falling edges of gate signal. Gate-in input must be considered a high impedance input and not a 50 Ω terminated input. In *Figure 2* is shown the block diagram of the Gate-In input. Gate-in input has a standard BNC connector. *Gate is active low: disable counting=TTL low, enable counting=TTL high.*




Figure 1: TTL OUT block diagram.




Figure 2 Gate-in Input block diagram

3.0 High count rate correction factor

It is well known that counting events with a detector that, after every detection, is blind for a fixed time, called dead-time (T_{deadtime}), causes an error in determining the true average counting rate. This happens because if n events are counted in a time window T_{measure} long, the correct observing time is not T_{measure} but T_{measure} reduced by the total blind time. The total blind-time can be easily calculated by multiplying the n counted events with T_{deadtime} . The actual count rate of the observed phenomena (CR_{actual}) can thus be calculated from the measured count rate (CR_{measure}) as follows:

$$CR_{\text{measured}} = \frac{n}{T_{\text{measure}}} \quad (\text{Eq. 1})$$

$$CR_{\text{actual}} = \frac{n}{T_{\text{measure}} - n \cdot T_{\text{dead-time}}} = \frac{CR_{\text{measured}}}{1 - CR_{\text{measured}} \cdot T_{\text{dead-time}}} \quad (\text{Eq. 2})$$

Although the correction is straightforward the closer the CR_{measured} is to the reciprocal of the dead-time, the higher is the correction and possibly the distortion to the measurement.

Dead-time correction can also be seen as a deviation from unity of the ratio between the actual incoming photon flux and the measured photon count rate. Given a dead-time, it is possible to calculate this ratio, called correction factor, as a function of the measured count rate (*Figure 3*).




Figure 3. Calculated correction factor for a PDM module with 62ns dead-time.

4.0 Optical Interface

SPAD sensor is placed in the centre of the front panel $\pm 0.25\text{mm}$ in both x and y directions. SPAD sensor surface is placed $2.6\text{mm} \pm 0.4\text{mm}$ from the window cap top glass surface.

